RadiPOP_API

Generated by Doxygen 1.9.2

1 Namespace Index 1
1.1 Packages
2 Hierarchical Index 3
2.1 Class Hierarchy
3 Class Index 5
3.1 Class List
4 File Index 7
4.1 File List
5 Namespace Documentation 9
5.1 config Namespace Reference
5.2 radipop_gui Namespace Reference
5.3 segmentation_utils Namespace Reference
5.3.1 Function Documentation
5.3.1.1 add_sobel_edges()
5.3.1.2 draw_region_outlines()
5.3.1.3 guess_bounds()
5.3.1.4 partition_at_threshold()
5.3.1.5 save_partition()
5.4 segmenter_flask_API Namespace Reference
5.4.1 Function Documentation
5.4.1.1 correctPartition()
5.4.1.2 dcm2png()
5.4.1.3 dcm2pngPreview()
5.4.1.4 drawOnMask()
5.4.1.5 extendLabels()
5.4.1.6 getMask()
5.4.1.7 highlightOrgan()
5.4.1.8 initialize()
5.4.1.9 labelOrgan()
5.4.1.10 postPickleGetMask()
5.4.1.11 saveMasks()
5.4.1.12 updateMask()
5.4.2 Variable Documentation
5.4.2.1 app
5.4.2.2 FLASK_PORT
5.4.2.3 FLAST_HOST
5.4.2.4 host
5.4.2.5 patients
5.4.2.6 port

6 Class Documentation	21
6.1 BaseConfig Class Reference	21
6.1.1 Member Data Documentation	21
6.1.1.1 SECRET_KEY	21
6.2 DevelopmentConfig Class Reference	21
6.2.1 Member Data Documentation	22
6.2.1.1 DEBUG	22
6.3 ProductionConfig Class Reference	22
6.4 RadiPopGUI Class Reference	22
6.4.1 Detailed Description	24
6.4.2 Constructor & Destructor Documentation	24
6.4.2.1init()	24
6.4.3 Member Function Documentation	25
6.4.3.1 clip_dcm()	25
6.4.3.2 color_mask()	25
6.4.3.3 correct_partition()	26
6.4.3.4 create_image_stream()	26
6.4.3.5 draw_on_image()	26
6.4.3.6 extend_labels()	27
6.4.3.7 extract_metadata_from_dcm()	27
6.4.3.8 find_organs()	27
6.4.3.9 highlight_np_label_array_to_png()	29
6.4.3.10 highlightOrgan()	30
6.4.3.11 labelMask()	30
6.4.3.12 np_label_array_to_png()	30
6.4.3.13 read_pickle_mask_to_np_label_array()	31
6.4.3.14 readPNG()	31
6.4.3.15 save_masks()	32
6.4.3.16 slice_dim()	32
6.4.3.17 update_mask_upon_slider_change()	32
6.4.3.18 writePillow2PNG()	33
6.4.4 Member Data Documentation	33
6.4.4.1 LABELS	33
6.4.4.2 last_clicked_x	33
6.4.4.3 last_clicked_y	33
6.4.4.4 LIVER_LABEL	34
6.4.4.5 masks	34
6.4.4.6 pathToSlices	34
6.4.4.7 patient_id	34
6.4.4.8 selected_pixel_value_of_label_mask	34
6.4.4.9 sliceCache	34
6.4.4.10 SPLEEN_LABEL	35

6.5 TestingConfig Class Reference	35
6.5.1 Member Data Documentation	35
6.5.1.1 DEBUG	35
7 File Documentation	37
7.1 /Users/lorenz/Desktop/temp/radipop/electron-react/segmenter_flask_API/segmenter_flask_API.py File Reference	37
7.2 /Users/lorenz/Desktop/temp/radipop/electron-react/segmenter_flask_API/utility/config.py File Reference	38
7.3 /Users/lorenz/Desktop/temp/radipop/electron-react/segmenter_flask_API/utility/radipop_gui.py File Reference	38
7.4 /Users/lorenz/Desktop/temp/radipop/electron-react/segmenter_flask_API/utility/segmentation_utils.py File Reference	38
Index	41

Namespace Index

1.1 Packages

Here are the packages with brief descriptions (if available):

config	9
radipop_gui	Ş
segmentation_utils	ξ
segmenter flask API	12

2 Namespace Index

Hierarchical Index

2.1 Class Hierarchy

This inheritance list is sorted roughly, but not completely, alphabetically:

BaseConfig	21
DevelopmentConfig	21
ProductionConfig	22
TestingConfig	35
RadiPopGUI	22

4 Hierarchical Index

Class Index

3.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:

BaseConfig	21
DevelopmentConfig	21
ProductionConfig	22
RadiPopGUI	
Bridge between the flask server/API and the RadiPOP segmenter:	22
TestingConfig	35

6 Class Index

File Index

4.1 File List

Here is a list of all files with brief descriptions:

/Users/lorenz/Desktop/temp/radipop/electron-react/segmenter_flask_API.py	37
/Users/lorenz/Desktop/temp/radipop/electron-react/segmenter_flask_API/utility/config.py	38
/Users/lorenz/Desktop/temp/radipop/electron-react/segmenter_flask_API/utility/radipop_gui.py	38
/Users/lorenz/Desktop/temp/radipop/electron-react/segmenter flask API/utility/segmentation utils.py	38

8 File Index

Namespace Documentation

5.1 config Namespace Reference

Classes

- class BaseConfig
- · class DevelopmentConfig
- · class ProductionConfig
- class TestingConfig

5.2 radipop_gui Namespace Reference

Classes

class RadiPopGUI

Bridge between the flask server/API and the RadiPOP segmenter:

5.3 segmentation_utils Namespace Reference

Functions

• def add_sobel_edges (mask, img)

Smooth edges Steps:

def draw_region_outlines (mask)

Color the mask light green.

def guess_bounds (regions_map, reference_map)

Guess the bounds/labels of the region based on reference region Guess the bounds/labels of the region based on reference region (generally neighboring slice).

• def partition_at_threshold (img, thresh, square_size, min_size, title=None, show_plot=True)

After some smoothing, calculate new mask for img Steps:

• def save_partition (mask, path)

5.3.1 Function Documentation

5.3.1.1 add_sobel_edges()

Smooth edges Steps:

- Edge filter image using the Canny algorithm.
- euclidean distance transform

Parameters

mask	mask corresponding to image
img	image corresponding to mask

Returns

mask with smoothed edges

5.3.1.2 draw_region_outlines()

```
\label{eq:continuous} \mbox{def segmentation\_utils.draw\_region\_outlines (} \\ \mbox{\it mask })
```

Color the mask light green.

Color the edges of the mask darker green.

Parameters

mask mask for which to color the outl	ines
---------------------------------------	------

Returns

mask with colored outlines

5.3.1.3 guess_bounds()

Guess the bounds/labels of the region based on reference region Guess the bounds/labels of the region based on reference region (generally neighboring slice).

Parameters

regions_map	mask to guess labels for
reference_map	reference mask (already labelled)

Returns

mask (labelled)

5.3.1.4 partition_at_threshold()

```
def segmentation_utils.partition_at_threshold (
    img,
    thresh,
    square_size,
    min_size,
    title = None,
    show_plot = True )
```

After some smoothing, calculate new mask for img Steps:

- · gaussian filter,
- · remove small objects,
- · greyscale morphological closing,
- euclidean distance transform

img	type numpy.ndarray
thres	Threshold value
min_size	Minimum size of an organ in the mask
squaresize	For greyscale morphological closing
title	Title of plot
show_plot	Show plot True/False

Returns

New binary mask (same size as img)

5.3.1.5 save partition()

5.4 segmenter_flask_API Namespace Reference

Functions

def correctPartition ()

Reveives index to slice/mask + coordinates--> returns partion corrected mask as PNG to client.

def dcm2png ()

Receive Paths to dcm files, converts them to PNG.

• def dcm2pngPreview ()

Receive Paths to dcm file, converts it to PNG.

def drawOnMask ()

Reveives index to slice/mask + x,y coordinates --> returns drawn on mask as PNG to client.

· def extendLabels ()

Reveives index to slice mask + label id --> returns highlighted mask as PNG to client.

def getMask ()

Reveives index to slice/mask --> returns mask stored on flask server as PNG to client.

• def highlightOrgan ()

Reveives index of slice + x,y coordinates --> returns highlighted mask as PNG to client.

• def initialize ()

Receive Paths to ordered slices, caches slices.

• def labelOrgan ()

Reveives index to slice mask + label id --> returns highlighted mask as PNG to client.

def postPickleGetMask ()

Receives path to pickle file --> returns mask as PNG to client.

• def saveMasks ()

Reveives path, saves all stored masks as pickle files to path --> returns output path.

def updateMask ()

Receives index of slice + slider values --> returns updated mask as PNG to client.

Variables

```
    app = Flask(__name___)
```

- int FLASK_PORT = 4041
- string FLAST_HOST = '0.0.0.0'
- host
- dictionary patients = {}

Dictionary which will hold for each patientID a RadiPopGUI object.

port

5.4.1 Function Documentation

5.4.1.1 correctPartition()

```
def segmenter_flask_API.correctPartition ( )
```

Reveives index to slice/mask + coordinates--> returns partion corrected mask as PNG to client.

Parameters

patientID	The ID of the patient
index	The index of the slice for which to mask should be updated
coordinates	array of coordinates of the form [x0,y0,x1,y1,,xn,yn]

Returns

JSON: {mask: byte stream} mask as transparent PNG as byte stream

Note: The coordinates array will be used to generate a line that cuts/divides the segmented organs.

Example handling of return image stream in js:

```
fetch(RadiPOP_states.FLASK_SERVER+"/labelOrgan", {
    method: 'POST',
    headers: { 'Content-Type': 'application/json'},
    body: JSON.stringify(data)
})
.then(function(response) { return response.json();})
.then(function(data) {
    bytestring = data["mask"];
    img = bytestring.split('\'')[1]
```

5.4.1.2 dcm2png()

```
def segmenter_flask_API.dcm2png ( )
```

Receive Paths to dcm files, converts them to PNG.

 $Included\ metadata\ IDs:\ "PatientID", "PatientBirthDate", "PatientName",\ "PatientAge", "PatientSex", "PatientName",\ "ContentDate"$

paths	An array with the paths to the slices
low_clip	lowest pixel value (Recommended:850)
high_clip	highest pixel value (Recommended: 1250)

Returns

JSON: {message: message, metadata: dictionary}

5.4.1.3 dcm2pngPreview()

```
def segmenter_flask_API.dcm2pngPreview ( )
```

Receive Paths to dcm file, converts it to PNG.

Included metadata IDs: "PatientID", "PatientBirthDate", "PatientName", "PatientAge", "PatientSex", "PatientName", " \leftarrow SliceThickness", "StudyID", "ContentDate"

Parameters

path	The paths to the dcm file
low_clip	lowest pixel value (Recommended:850)
high_clip	highest pixel value (Recommended: 1250)

Returns

JSON: {slice: slice stream} slice as PNG as byte stream

5.4.1.4 drawOnMask()

```
def segmenter_flask_API.drawOnMask ( )
```

Reveives index to slice/mask + x,y coordinates --> returns drawn on mask as PNG to client.

Parameters

patientID	The ID of the patient
index	The index of the slice for which to mask should be updated
coordinates	array of coordinates of the form [x0,y0,x1,y1,,xn,yn]

Returns

JSON: {mask: byte stream} mask as transparent PNG as byte stream

Note: The coordinates array will be used to draw a line on the mask.

Example handling of return image stream in js:

```
fetch(RadiPOP_states.FLASK_SERVER+"/labelOrgan", {
    method: 'POST',
```

```
headers: { 'Content-Type': 'application/json'},
   body: JSON.stringify(data)
})
.then(function(response) { return response.json();})
.then(function(data) {
bytestring = data["mask"];
img = bytestring.split('\'')[1]
```

5.4.1.5 extendLabels()

```
def segmenter_flask_API.extendLabels ( )
```

Reveives index to slice mask + label id --> returns highlighted mask as PNG to client.

Parameters

patientID	The ID of the patient
index	The index of the slice for which to mask should be updated
left	Extend labeling up to index-left
right	Extend labeling up to index+right

Returns

```
JSON: {left_most_idx: idx, right_most_idx: idx}
```

Note: The left_most_idx and right_most_idx correspond to the indices of the slices up to which the labeling has been extended. After the the function has finished use the function API's function /getMask to update the masks in your GUI. Example in js:

```
for (let index=parseInt(data["left_most_idx"]); index<parseInt(data["right_most_idx"])+1; index++) {
    $.post(FLASK_SERVER+"/getMask", {
        javascript_data: JSON.stringify({patienID: id, index: idx})
    })
}</pre>
```

5.4.1.6 getMask()

```
def segmenter_flask_API.getMask ( )
```

Reveives index to slice/mask --> returns mask stored on flask server as PNG to client.

patientID	The ID of the patient
index	The index of the slice for which to mask should be updated

Returns

JSON: {mask: byte stream} mask as transparent PNG as byte stream

Example handling of return image stream in js:

```
fetch(RadiPOP_states.FLASK_SERVER+"/labelOrgan", {
    method: 'POST',
    headers: { 'Content-Type': 'application/json'},
    body: JSON.stringify(data)
})
.then(function(response) { return response.json();})
.then(function(data) {
    bytestring = data["mask"];
    img = bytestring.split('\'')[1]
```

5.4.1.7 highlightOrgan()

```
def segmenter_flask_API.highlightOrgan ( )
```

Reveives index of slice + x,y coordinates --> returns highlighted mask as PNG to client.

Parameters

patientID	The ID of the patient
index	The index of the slice for which to mask should be updated
X	relative x coordinates (0<=x<=1)
у	relative y coordinates (0<=y<=1)

Returns

JSON: {mask: byte stream} mask as transparent PNG as byte stream

Example handling of return image stream in js:

```
fetch(RadiPOP_states.FLASK_SERVER+"/labelOrgan", {
    method: 'POST',
    headers: { 'Content-Type': 'application/json'},
    body: JSON.stringify(data)
})
.then(function(response) { return response.json();})
.then(function(data) {
    bytestring = data["mask"];
    img = bytestring.split('\'')[1]
```

5.4.1.8 initialize()

```
def segmenter_flask_API.initialize ( )
```

Receive Paths to ordered slices, caches slices.

Parameters

patientID	The ID of the patient. Must be unique!
paths	An array with the paths to the slices

Returns

JSON: {message: message}

Note: Paths to slices !!!MUST BE ORDERED!!! 0,1,..,n

5.4.1.9 labelOrgan()

```
def segmenter_flask_API.labelOrgan ( )
```

Reveives index to slice mask + label id --> returns highlighted mask as PNG to client.

Parameters

patientID	The ID of the patient
index	The index of the slice for which to mask should be updated
label	Label of organ (1 for liver, 2 for spleen, 0 nothing, >2 other organ)

Returns

JSON: {mask: byte stream} mask as transparent PNG as byte stream

Example handling of return image stream in js:

```
fetch(RadiPOP_states.FLASK_SERVER+"/labelOrgan", {
    method: 'POST',
    headers: { 'Content-Type': 'application/json'},
    body: JSON.stringify(data)
})
.then(function(response) { return response.json();})
.then(function(data) {
    bytestring = data["mask"];
    img = bytestring.split('\'')[1]
```

5.4.1.10 postPickleGetMask()

```
def segmenter_flask_API.postPickleGetMask ( )
```

Receives path to pickle file --> returns mask as PNG to client.

patientID	The ID of the patient
index	The index of the slice the mask refers to
Generaled by Do	xyJee path to the mask file

Returns

JSON: {mask: byte stream} mask as transparent PNG as byte stream

Example handling of return image stream in js:

```
fetch(RadiPOP_states.FLASK_SERVER+"/labelOrgan", {
    method: 'POST',
    headers: { 'Content-Type': 'application/json'},
    body: JSON.stringify(data)
})
.then(function(response) { return response.json();})
.then(function(data) {
    bytestring = data["mask"];
    img = bytestring.split('\'')[1]
```

5.4.1.11 saveMasks()

```
def segmenter_flask_API.saveMasks ( )
```

Reveives path, saves all stored masks as pickle files to path --> returns output path.

Parameters

patientID	The ID of the patient
index	The index of the slice for which to mask should be updated

Returns

JSON: {outdir: path} path/directory to which the pickle files were written

5.4.1.12 updateMask()

```
def segmenter_flask_API.updateMask ( )
```

Receives index of slice + slider values --> returns updated mask as PNG to client.

patientID	The ID of the patient
index	The index of the slice for which to mask should be updated
liver-intensity-slider	Slider value for liver intesity
bone-intensity-slider	Slider value for bone intesity
blood-vessel-intensity-slider	Slider value for blood-vessel intesity

Returns

JSON: {mask: byte stream} mask as transparent PNG as byte stream

Example handling of return image stream in js:

```
fetch(RadiPOP_states.FLASK_SERVER+"/labelOrgan", {
    method: 'POST',
    headers: { 'Content-Type': 'application/json'},
    body: JSON.stringify(data)
})
.then(function(response) { return response.json();})
.then(function(data) {
    bytestring = data["mask"];
    img = bytestring.split('\'')[1]
```

5.4.2 Variable Documentation

5.4.2.1 app

```
app = Flask(__name___)
```

5.4.2.2 FLASK_PORT

```
int FLASK_PORT = 4041
```

5.4.2.3 FLAST_HOST

```
FLAST_HOST = '0.0.0.0'
```

5.4.2.4 host

host

5.4.2.5 patients

```
dictionary patients = {}
```

Dictionary which will hold for each patientID a RadiPopGUI object.

Patients are added by the API's /initialize function


5.4.2.6 port

port

Class Documentation

6.1 BaseConfig Class Reference

Inheritance diagram for BaseConfig:



Static Public Attributes

• SECRET_KEY = os.getenv('SECRET_KEY', 'REPLACE ME')

6.1.1 Member Data Documentation

6.1.1.1 SECRET_KEY


```
SECRET_KEY = os.getenv('SECRET_KEY', 'REPLACE ME') [static]
```

The documentation for this class was generated from the following file:

 $\bullet \ / Users/lorenz/Desktop/temp/radipop/electron-react/segmenter_flask_API/utility/config.py$

6.2 DevelopmentConfig Class Reference

Inheritance diagram for DevelopmentConfig:



22 Class Documentation

Static Public Attributes

• bool DEBUG = True

6.2.1 Member Data Documentation

6.2.1.1 DEBUG


```
bool DEBUG = True [static]
```

The documentation for this class was generated from the following file:

/Users/lorenz/Desktop/temp/radipop/electron-react/segmenter_flask_API/utility/config.py

6.3 ProductionConfig Class Reference

Inheritance diagram for ProductionConfig:



Additional Inherited Members

The documentation for this class was generated from the following file:

/Users/lorenz/Desktop/temp/radipop/electron-react/segmenter_flask_API/utility/config.py

6.4 RadiPopGUI Class Reference

Bridge between the flask server/API and the RadiPOP segmenter:

Public Member Functions

def __init__ (self, patient_id)

Class constructor:

def extend_labels (self, cur_idx, left_extend, right_extend)

Extend labels from current slice to neighbouring slices.

def highlight_np_label_array_to_png (self, mask_np_array, highlight)

Takes an numpy label array dim(n,m,1) and returns a RGBA pillow image dim(n,m,4)

def highlightOrgan (self, slice_idx, x, y)

Highlights regions of the mask (organs) that were clicked on by user.

def labelMask (self, slice_idx, newlabel)

Labels mask at given index at previously selected region.

• def save_masks (self, path)

Saves masks as pickle file to given path.

• def slice dim (self, index)

Returns dimensions of slice images (x,y)

Static Public Member Functions

def clip_dcm (dcm_file, clip_low=850, clip_high=1250)

Read dicom image (.dcm), clips it and returns it as a grey scale PNG.

• def color mask (mask np array)

Takes an numpy label array dim(n,m,1) and returns a color mask dim(n,m,4)

def correct partition (image)

Convert PNG mask to 1 channelled label mask.

• def create_image_stream (img)

Returns base64 bytestream for given input image.

def draw_on_image (coordinates, img, correctionMode=False)

Draws a point or lines on given PNG image.

def extract_metadata_from_dcm (dcm_file)

Read dicom image (.dcm) and extract metadata Extracted metadata IDs: "PatientID", "PatientBirthDate", "Patient ← Name", "PatientAge", "PatientSex", "PatientName", "SliceThickness", "StudyID", "ContentDate".

• def find_organs (img, bones_thresh, blood_vessels_thresh, liver_thresh)

Uses three threshold values to find organs.

def np_label_array_to_png (mask_np_array)

Takes an numpy label array dim(n,m,1) and returns a RGBA pillow image dim(n,m,4)

· def read pickle mask to np label array (path)

Opens mask pickle file and returns it as a np.array.

• def readPNG (path)

Reads an image (e.g.

def update_mask_upon_slider_change (image, bone_intensity, blood_vessel_intensity, liver_intensity)

Sets threshold for liver intensity.

• def writePillow2PNG (img, outfile)

24 Class Documentation

Public Attributes

· last clicked x

The x-coordinate of the region that was last selected/clicked on.

last_clicked_y

The y-coordinate of the region that was last selected/clicked on.

· masks

Dictionary: key: mask index, value: mask numpy array dim(n,m,1)

· pathToSlices

List contaning the path to the png slice files

· patient id

ID of the patient.

• selected_pixel_value_of_label_mask

The label of the region that was last selected/clicked on.

sliceCache

Dictionary: key: slice index, value: slice PNG.

Static Public Attributes

```
• list LABELS = [SPLEEN_LABEL,SPLEEN_LABEL]
```

list with labels

• int LIVER LABEL = 1

Regions in mask with this label value are considered liver.

• int SPLEEN_LABEL = 2

Regions in mask with this label value are considered spleen.

6.4.1 Detailed Description

Bridge between the flask server/API and the RadiPOP segmenter:

Note:

- For each patient one object of this class should be instantiated.
- · An object of this class contains all the slices and masks associated with the patient
- · This class also contains static utility functions
- This class is the bridge between the flask server/API and the RadiPOP segmenter (segmentations_utils)

6.4.2 Constructor & Destructor Documentation

Class constructor:

Parameters

patient←	The ID of the patient. Must be unique!
_id	

Note:

- · For each patient one object of this class should be instantiated.
- · An object of this class contains all the slices and masks associated with the patient

6.4.3 Member Function Documentation

6.4.3.1 clip_dcm()

Read dicom image (.dcm), clips it and returns it as a grey scale PNG.

```
@param dcm_file Path to .dcm file
@param clip_low lowest pixel value
@param clip_high highest pixel value
@return tuple(L (grey scale) Pillow Image, slice index)
```

6.4.3.2 color_mask()

Takes an numpy label array dim(n,m,1) and returns a color mask dim(n,m,4)

Parameters

```
mask_np_array | numpy label array dim(n,m,1)
```

Returns

color mask image dim(n,m,4) --> RGBA

```
Default liver color is red (LIVER_LABEL)Default spleen color is blue (SPLEEN_LABEL)Default for other regions is green
```

26 Class Documentation

6.4.3.3 correct_partition()

Convert PNG mask to 1 channelled label mask.

Parameters

```
image | Image (e.g.: RGBA PNG pillow)
```

Returns

label mask

6.4.3.4 create_image_stream()

```
\begin{tabular}{ll} def & create_image_stream & ( & \\ & img & ) & [static] \\ \end{tabular}
```

Returns base64 bytestream for given input image.

Parameters

```
img Image (e.g. Pillow PNG)
```

Returns

img_base64 stream

6.4.3.5 draw_on_image()

Draws a point or lines on given PNG image.

coordinates	List of the form [x0,y0,x1,y1,,xn,yn]
img	Image (e.g.: RGBA PNG pillow)
correctionMode	True/False If true the drawn line acts as an eraser (dividing organs). If false acentared புறைப்படி drawn on the image. DEFAULT: False

The modifications are made directly on the provided image. No return value

6.4.3.6 extend labels()

Extend labels from current slice to neighbouring slices.

Parameters

cur_idx	Index of reference slice
left_extend	Number of slices to extend the labeling to below reference slice
right_extend	Number of slices to extend the labeling to above reference slice

Returns

(left_most_idx,right_most_idx) The index of the outermost slices to which the labeling was extended

Extends labels left and right from current slice How far the labels are extended is taken from left and right expansion bounds

6.4.3.7 extract_metadata_from_dcm()

Read dicom image (.dcm) and extract metadata Extracted metadata IDs: "PatientID", "PatientBirthDate", "Patient Wame", "PatientAge", "PatientSex", "PatientName", "SliceThickness", "StudyID", "ContentDate".

Parameters

```
dcm_file Path to .dcm file
```

Returns

dictionary with metadata information

6.4.3.8 find_organs()

```
def find_organs (
          img,
          bones_thresh,
```

28 Class Documentation

blood_vessels_thresh,
liver_thresh) [static]

Uses three threshold values to find organs.

Parameters

image	Image (e.g.: RGBA PNG pillow)
bones_thresh	bones threshold: [threshold, square_size , min_size]
blood_vessels_thresh	blood vessels threshold: [threshold, square_size , min_size]
liver_thresh	liver threshold: [threshold, square_size, min_size]

Returns

New labelled mask (same size as slice)

```
The algorithm is:

- After some smoothing, remove every pixel above bones threshold from the image.

- After some smoothing, remove every pixel above blood vessel threshold.

- Everything that then remains above liver threshold is called an organ.

- Use contiguous area divisions to roughly split into organs.
```

6.4.3.9 highlight_np_label_array_to_png()

Takes an numpy label array dim(n,m,1) and returns a RGBA pillow image dim(n,m,4)

Parameters

mask_np_array	numpy label array dim(n,m,1) @highlight Highlight regions where highlight==label in brighter	
	color	

Returns

RGBA pillow image dim(n,m,4)

Turns a labelled mask into a transparent (RGBA) PNG.

- Default liver color is red (LIVER_LABEL)
- Default spleen color is blue (SPLEEN_LABEL)
- · Default for other regions is green

30 Class Documentation

6.4.3.10 highlightOrgan()

```
def highlightOrgan (
    self,
    slice_idx,
    x,
    y )
```

Highlights regions of the mask (organs) that were clicked on by user.

Parameters

slice_idx	Index of mask/slice to be highlighted
X	x-coordinate of slice (in pixels)
У	y-coordinate of slice (in pixels)

Returns

tuple(Mask where the region specified by x and y is highlighted in brigther colors, pixel_value)

6.4.3.11 labelMask()

Labels mask at given index at previously selected region.

Parameters

slice_idx	Index of mask/slice to be labelled
label	Label to be assigned to previously selected region (either LIVER_LABEL or SPLEEN_LABEL or -1
	(for removeing label))

Returns

Mask with new label

Note: It is expected that the client has before highlighted an organ with the function self.highlightOrgan(). This determines the region/organ that will be labelled with label.

6.4.3.12 np_label_array_to_png()

Takes an numpy label array dim(n,m,1) and returns a RGBA pillow image dim(n,m,4)

Parameters

mask_np_array	numpy label array dim(n,m,1) @highlight Highlight regions where highlight==label in brighter	
	color	

Returns

RGBA pillow image dim(n,m,4)

Turns a labelled mask into a transparent (RGBA) PNG.

- Default liver color is red (LIVER_LABEL)
- Default spleen color is blue (SPLEEN_LABEL)
- Default for other regions is green

6.4.3.13 read_pickle_mask_to_np_label_array()

```
\label{local_def} \mbox{def read_pickle_mask_to_np_label_array (} \\ path \ ) \ \ [static]
```

Opens mask pickle file and returns it as a np.array.

Parameters

path	Path to pickle file

Returns

numpy array of mask

6.4.3.14 readPNG()

```
\label{eq:continuous_path} \begin{array}{ll} \text{def readPNG (} \\ & path \text{ )} & [\text{static}] \end{array}
```

Reads an image (e.g.

: PNG file) to numpy array

Parameters

path Path to image

32 Class Documentation

Returns

numpy array of image

6.4.3.15 save_masks()

```
def save_masks (
          self,
          path )
```

Saves masks as pickle file to given path.

Parameters

Masks are written as .p (pickle) files

6.4.3.16 slice_dim()

```
def slice_dim (
          self,
          index )
```

Returns dimensions of slice images (x,y)

Returns

(x,y) Dimensions of slices

6.4.3.17 update_mask_upon_slider_change()

Sets threshold for liver intensity.

Parameters

image	Image (e.g.: RGBA PNG pillow)	
bones_thresh	bones threshold	
blood_vessels_thresh	blood vessels threshold	
liver_thresh	liver threshold	

Returns

New labelled mask

Steps:

- · Sets thresholds for current slice
- Runs self.find_organs on current slice with new thresholds

6.4.3.18 writePillow2PNG()

6.4.4 Member Data Documentation

6.4.4.1 LABELS

```
list LABELS = [SPLEEN_LABEL, SPLEEN_LABEL] [static]
```

list with labels

6.4.4.2 last_clicked_x

```
last_clicked_x
```

The x-coordinate of the region that was last selected/clicked on.

6.4.4.3 last_clicked_y

```
last_clicked_y
```

The y-coordinate of the region that was last selected/clicked on.

34 Class Documentation

6.4.4.4 LIVER_LABEL

```
int LIVER_LABEL = 1 [static]
```

Regions in mask with this label value are considered liver.

6.4.4.5 masks

masks

Dictionary: key: mask index, value: mask numpy array dim(n,m,1)

6.4.4.6 pathToSlices

pathToSlices

List contaning the path to the png slice files

6.4.4.7 patient_id

patient_id

ID of the patient.

6.4.4.8 selected_pixel_value_of_label_mask

```
selected_pixel_value_of_label_mask
```

The label of the region that was last selected/clicked on.

6.4.4.9 sliceCache

sliceCache

Dictionary: key: slice index, value: slice PNG.

6.4.4.10 SPLEEN_LABEL

```
int SPLEEN_LABEL = 2 [static]
```


Regions in mask with this label value are considered spleen.

The documentation for this class was generated from the following file:

• /Users/lorenz/Desktop/temp/radipop/electron-react/segmenter_flask_API/utility/radipop_gui.py

6.5 TestingConfig Class Reference

Inheritance diagram for TestingConfig:



Static Public Attributes

• bool DEBUG = True

6.5.1 Member Data Documentation

6.5.1.1 DEBUG

```
bool DEBUG = True [static]
```

The documentation for this class was generated from the following file:

• /Users/lorenz/Desktop/temp/radipop/electron-react/segmenter_flask_API/utility/config.py

36 Class Documentation

Chapter 7

File Documentation

7.1 /Users/lorenz/Desktop/temp/radipop/electron-react/segmenter_← flask_API/segmenter_flask_API.py File Reference

Namespaces

· namespace segmenter flask API

Functions

def correctPartition ()

Reveives index to slice/mask + coordinates--> returns partion corrected mask as PNG to client.

• def dcm2png ()

Receive Paths to dcm files, converts them to PNG.

def dcm2pngPreview ()

Receive Paths to dcm file, converts it to PNG.

def drawOnMask ()

Reveives index to slice/mask + x,y coordinates --> returns drawn on mask as PNG to client.

• def extendLabels ()

Reveives index to slice mask + label id --> returns highlighted mask as PNG to client.

· def getMask ()

Reveives index to slice/mask --> returns mask stored on flask server as PNG to client.

• def highlightOrgan ()

Reveives index of slice + x,y coordinates --> returns highlighted mask as PNG to client.

• def initialize ()

Receive Paths to ordered slices, caches slices.

• def labelOrgan ()

Reveives index to slice mask + label id --> returns highlighted mask as PNG to client.

def postPickleGetMask ()

Receives path to pickle file --> returns mask as PNG to client.

• def saveMasks ()

Reveives path, saves all stored masks as pickle files to path --> returns output path.

• def updateMask ()

Receives index of slice + slider values --> returns updated mask as PNG to client.

38 File Documentation

Variables

```
app = Flask(__name__)
int FLASK_PORT = 4041
string FLAST_HOST = '0.0.0.0'
host
dictionary patients = {}

Dictionary which will hold for each patientID a RadiPopGUI object.
port
```

7.2 /Users/lorenz/Desktop/temp/radipop/electron-react/segmenter_← flask API/utility/config.py File Reference

Classes

- · class BaseConfig
- · class DevelopmentConfig
- class ProductionConfig
- · class TestingConfig

Namespaces

· namespace config

7.3 /Users/lorenz/Desktop/temp/radipop/electron-react/segmenter_← flask_API/utility/radipop_gui.py File Reference

Classes

class RadiPopGUI

Bridge between the flask server/API and the RadiPOP segmenter:

Namespaces

- namespace radipop_gui
- 7.4 /Users/lorenz/Desktop/temp/radipop/electron-react/segmenter_← flask_API/utility/segmentation_utils.py File Reference

Namespaces

• namespace segmentation_utils

Functions

• def add_sobel_edges (mask, img)

Smooth edges Steps:

• def draw_region_outlines (mask)

Color the mask light green.

• def guess_bounds (regions_map, reference_map)

Guess the bounds/labels of the region based on reference region Guess the bounds/labels of the region based on reference region (generally neighboring slice).

• def partition_at_threshold (img, thresh, square_size, min_size, title=None, show_plot=True)

After some smoothing, calculate new mask for img Steps:

def save_partition (mask, path)

40 File Documentation

Index

```
/Users/lorenz/Desktop/temp/radipop/electron-react/segmentotrafitasknetatadistaginemitedcflask_API.py,
                                                        RadiPopGUI, 27
/Users/lorenz/Desktop/temp/radipop/electron-react/segmenter flask API/utility/config.py,
                                                   find_organs
/Users/lorenz/Desktop/temp/radipop/electron-react/segmenter flastipaps/lity/radipop gui.py,
                                                   FLASK PORT
FLAST_HOST
                                                       segmenter_flask_API, 19
 init
    RadiPopGUI, 24
                                                   getMask
                                                       segmenter_flask_API, 15
add sobel edges
                                                   guess_bounds
    segmentation utils, 10
                                                       segmentation_utils, 10
app
    segmenter flask API, 19
                                                   highlight_np_label_array_to_png
                                                       RadiPopGUI, 29
BaseConfig, 21
                                                   highlightOrgan
    SECRET KEY, 21
                                                       RadiPopGUI, 29
clip dcm
                                                       segmenter flask API, 16
    RadiPopGUI, 25
                                                   host
color mask
                                                       segmenter_flask_API, 19
    RadiPopGUI, 25
                                                   initialize
config, 9
                                                       segmenter_flask_API, 16
correct_partition
    RadiPopGUI, 26
                                                   labelMask
correctPartition
                                                       RadiPopGUI, 30
    segmenter flask API, 13
                                                   labelOrgan
create image stream
                                                       segmenter_flask_API, 17
    RadiPopGUI, 26
                                                   LABELS
dcm2png
                                                        RadiPopGUI, 33
    segmenter flask API, 13
                                                   last_clicked_x
dcm2pngPreview
                                                       RadiPopGUI, 33
    segmenter_flask_API, 14
                                                   last_clicked_y
DEBUG
                                                        RadiPopGUI, 33
    DevelopmentConfig, 22
                                                   LIVER LABEL
    TestingConfig, 35
                                                       RadiPopGUI, 33
DevelopmentConfig, 21
                                                   masks
    DEBUG, 22
                                                       RadiPopGUI, 34
draw on image
    RadiPopGUI, 26
                                                   np_label_array_to_png
draw region outlines
                                                       RadiPopGUI, 30
    segmentation utils, 10
drawOnMask
                                                   partition_at_threshold
    segmenter_flask_API, 14
                                                       segmentation_utils, 11
                                                   pathToSlices
extend labels
                                                       RadiPopGUI, 34
    RadiPopGUI, 27
                                                   patient id
extendLabels
                                                       RadiPopGUI, 34
    segmenter_flask_API, 15
```

42 INDEX

patients	app, 19
segmenter_flask_API, 19	correctPartition, 13
port	dcm2png, 13
segmenter_flask_API, 19	dcm2pngPreview, 14
postPickleGetMask	drawOnMask, 14
segmenter_flask_API, 17	extendLabels, 15
ProductionConfig, 22	FLASK_PORT, 19
	FLAST_HOST, 19
radipop_gui, 9	getMask, 15
RadiPopGUI, 22	highlightOrgan, 16
init, 24	host, 19
clip_dcm, 25	initialize, 16
color_mask, 25	labelOrgan, 17
correct_partition, 26	patients, 19
create_image_stream, 26	port, 19
draw_on_image, 26	postPickleGetMask, 17
extend_labels, 27	saveMasks, 18
extract_metadata_from_dcm, 27	updateMask, 18
find_organs, 27	selected_pixel_value_of_label_mask
highlight_np_label_array_to_png, 29	RadiPopGUI, 34
highlightOrgan, 29	slice_dim
labelMask, 30	RadiPopGUI, 32
LABELS, 33	sliceCache
last_clicked_x, 33	RadiPopGUI, 34
last_clicked_y, 33	SPLEEN_LABEL
LIVER_LABEL, 33	RadiPopGUI, 34
masks, 34	
np_label_array_to_png, 30	TestingConfig, 35
pathToSlices, 34	DEBUG, 35
patient_id, 34	
read_pickle_mask_to_np_label_array, 31	update_mask_upon_slider_change
readPNG, 31	RadiPopGUI, 32
save_masks, 32	updateMask segmenter flask API, 18
selected_pixel_value_of_label_mask, 34	segmenter_nask_AFI, To
slice_dim, 32	writePillow2PNG
sliceCache, 34	RadiPopGUI, 33
SPLEEN_LABEL, 34	
update_mask_upon_slider_change, 32 writePillow2PNG, 33	
read_pickle_mask_to_np_label_array	
RadiPopGUI, 31	
readPNG	
RadiPopGUI, 31	
agya magka	
save_masks	
RadiPopGUI, 32	
save_partition	
segmentation_utils, 12 saveMasks	
segmenter_flask_API, 18	
SECRET_KEY Reso Config. 21	
BaseConfig, 21	
segmentation_utils, 9	
add_sobel_edges, 10	
draw_region_outlines, 10	
guess_bounds, 10	
partition_at_threshold, 11	
save_partition, 12	
segmenter_flask_API, 12	