

RadiPOP_API

Generated by Doxygen 1.9.2

1 RadiPOP segmenter backend	1
1.1 Description	1
2 Namespace Index	3
2.1 Packages	3
3 Class Index	5
3.1 Class List	5
4 File Index	7
4.1 File List	7
5 Namespace Documentation	9
5.1 radipop_gui Namespace Reference	9
5.1.1 Detailed Description	9
5.1.2 Description	9
5.1.3 Author(s)	9
5.2 segmentation_utils Namespace Reference	10
5.2.1 Detailed Description	10
5.2.2 Author(s)	10
5.2.3 Function Documentation	10
5.2.3.1 add_sobel_edges()	10
5.2.3.2 draw_region_outlines()	11
5.2.3.3 guess_bounds()	11
5.2.3.4 partition_at_threshold()	11
5.2.3.5 save_partition()	12
5.3 segmenter_flask_API Namespace Reference	12
5.3.1 Detailed Description	13
5.3.2 Description	13
5.3.3 Author(s)	13
5.3.4 Function Documentation	14
5.3.4.1 correctPartition()	14
5.3.4.2 dcm2png()	14
5.3.4.3 dcm2pngPreview()	15
5.3.4.4 drawOnMask()	15
5.3.4.5 extendLabels()	16
5.3.4.6 getMask()	16
5.3.4.7 highlightOrgan()	17
5.3.4.8 initialize()	17
5.3.4.9 labelOrgan()	18
5.3.4.10 postPickleGetMask()	18
5.3.4.11 saveMasks()	19
5.3.4.12 updateMask()	19
5.3.5 Variable Documentation	20

5.3.5.1 app	20
5.3.5.2 FLASK_PORT	20
5.3.5.3 FLAST_HOST	20
5.3.5.4 host	20
5.3.5.5 patients	20
5.3.5.6 port	20
6 Class Documentation	21
6.1 RadiPopGUI Class Reference	21
6.1.1 Detailed Description	22
6.1.2 Constructor & Destructor Documentation	23
6.1.2.1 __init__()	23
6.1.3 Member Function Documentation	23
6.1.3.1 clip_dcm()	23
6.1.3.2 color_mask()	23
6.1.3.3 correct_partition()	24
6.1.3.4 create_image_stream()	24
6.1.3.5 draw_on_image()	25
6.1.3.6 extend_labels()	25
6.1.3.7 extract_metadata_from_dcm()	25
6.1.3.8 find_organs()	26
6.1.3.9 highlight_np_label_array_to_png()	26
6.1.3.10 highlightOrgan()	27
6.1.3.11 labelMask()	27
6.1.3.12 np_label_array_to_png()	28
6.1.3.13 read_pickle_mask_to_np_label_array()	28
6.1.3.14 readPNG()	29
6.1.3.15 save_masks()	29
6.1.3.16 slice_dim()	29
6.1.3.17 update_mask_upon_slider_change()	30
6.1.3.18 writePillow2PNG()	30
6.1.4 Member Data Documentation	30
6.1.4.1 LABELS	30
6.1.4.2 last_clicked_x	31
6.1.4.3 last_clicked_y	31
6.1.4.4 LIVER_LABEL	31
6.1.4.5 masks	31
6.1.4.6 pathToSlices	31
6.1.4.7 patient_id	31
6.1.4.8 selected_pixel_value_of_label_mask	32
6.1.4.9 sliceCache	32
6.1.4.10 SPLEEN_LABEL	32

7 File Documentation	33
7.1 /Users/lorenz/Desktop/radipop/electron-react/segmenter_flask_API/segmenter_flask_API.py File Reference	33
7.2 /Users/lorenz/Desktop/radipop/electron-react/segmenter_flask_API/utility/radipop_gui.py File Reference	34
7.3 /Users/lorenz/Desktop/radipop/electron-react/segmenter_flask_API/utility/segmentation_utils.py File Reference	34
Index	35

Chapter 1

RadiPOP segmenter backend

1.1 Description

This python code constitutes the backend of the RadiPOP segmenter.

- Contains code for handling communication with frontend - Flask server ([segmenter_flask_API.py](#))
- Contains python class for storing patient information, handling flask requests and calculations ([utility\radipop_gui.py](#))
- Contains code from the original segmenter: Segmentation algorithm, applying thresholds ... ([utility\segmentation_utils.py](#))

Chapter 2

Namespace Index

2.1 Packages

Here are the packages with brief descriptions (if available):

radipop_gui	Contains python class for storing patient information, handling flask requests and calculations	9
segmentation_utils	Contains code from the original segmenter: Segmentation algorithm, applying thresholds	10
segmenter_flask_API	Flask server for RadiPOP segmenter: Bridge between frontend and backend	12

Chapter 3

Class Index

3.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:

[RadiPopGUI](#)

Bridge between the flask server/API and the RadiPOP segmenter: [21](#)

Chapter 4

File Index

4.1 File List

Here is a list of all files with brief descriptions:

/Users/lorenz/Desktop/radipop/electron-react/segmenter_flask_API/ segmenter_flask_API.py	33
/Users/lorenz/Desktop/radipop/electron-react/segmenter_flask_API/utility/ radipop_gui.py	34
/Users/lorenz/Desktop/radipop/electron-react/segmenter_flask_API/utility/ segmentation_utils.py	34

Chapter 5

Namespace Documentation

5.1 radipop_gui Namespace Reference

Contains python class for storing patient information, handling flask requests and calculations.

Classes

- class [RadiPopGUI](#)

Bridge between the flask server/API and the RadiPOP segmenter:

5.1.1 Detailed Description

Contains python class for storing patient information, handling flask requests and calculations.

5.1.2 Description

Bridge between the flask server/API and the RadiPOP segmenter:

- For each patient one object of this class should be instantiated.
- An object of this class contains all the slices and masks associated with the patient
- This class also contains static utility functions
- This class is the bridge between the flask server/API and the RadiPOP segmenter (segmentations_utils)

5.1.3 Author(s)

- Created by Lorenz Perschy (2021, 2022)

5.2 segmentation_utils Namespace Reference

Contains code from the original segmenter: Segmentation algorithm, applying thresholds ...

Functions

- def [add_sobel_edges](#) (mask, img)
Smooth edges Steps:
- def [draw_region_outlines](#) (mask)
Color the mask light green.
- def [guess_bounds](#) (regions_map, reference_map)
Guess the bounds/labels of the region based on reference region Guess the bounds/labels of the region based on reference region (generally neighboring slice).
- def [partition_at_threshold](#) (img, thresh, square_size, min_size, title=None, show_plot=True)
After some smoothing, calculate new mask for img Steps:
- def [save_partition](#) (mask, path)

5.2.1 Detailed Description

Contains code from the original segmenter: Segmentation algorithm, applying thresholds ...

5.2.2 Author(s)

Menche Lab

5.2.3 Function Documentation

5.2.3.1 [add_sobel_edges\(\)](#)

```
def segmentation_utils.add_sobel_edges (
    mask,
    img )
```

Smooth edges Steps:

- Edge filter image using the Canny algorithm.
- euclidean distance transform

Parameters

<i>mask</i>	mask corresponding to image
<i>img</i>	image corresponding to mask

Returns

mask with smoothed edges

5.2.3.2 draw_region_outlines()

```
def segmentation_utils.draw_region_outlines (
    mask )
```

Color the mask light green.

Color the edges of the mask darker green.

Parameters

<i>mask</i>	mask for which to color the outlines
-------------	--------------------------------------

Returns

mask with colored outlines

5.2.3.3 guess_bounds()

```
def segmentation_utils.guess_bounds (
    regions_map,
    reference_map )
```

Guess the bounds/labels of the region based on reference region Guess the bounds/labels of the region based on reference region (generally neighboring slice).

Parameters

<i>regions_map</i>	mask to guess labels for
<i>reference_map</i>	reference mask (already labelled)

Returns

mask (labelled)

5.2.3.4 partition_at_threshold()

```
def segmentation_utils.partition_at_threshold (
    img,
```

```

    thresh,
    square_size,
    min_size,
    title = None,
    show_plot = True )

```

After some smoothing, calculate new mask for img Steps:

- gaussian filter,
- remove small objects,
- greyscale morphological closing,
- euclidean distance transform

Parameters

<i>img</i>	type numpy.ndarray
<i>thres</i>	Threshold value
<i>min_size</i>	Minimum size of an organ in the mask
<i>squaresize</i>	For greyscale morphological closing
<i>title</i>	Title of plot
<i>show_plot</i>	Show plot True/False

Returns

New binary mask (same size as img)

5.2.3.5 save_partition()

```

def segmentation_utils.save_partition (
    mask,
    path )

```

5.3 segmenter_flask_API Namespace Reference

Flask server for RadiPOP segmenter: Bridge between frontend and backend.

Functions

- def [correctPartition](#) ()
Reveives index to slice/mask + coordinates--> returns partion corrected mask as PNG to client.
- def [dcm2png](#) ()
Receive Paths to dcm files, converts them to PNG.
- def [dcm2pngPreview](#) ()

- Receive Paths to dcm file, converts it to PNG.*
 - def `drawOnMask` ()
 - Reveives index to slice/mask + x,y coordinates --> returns drawn on mask as PNG to client.*
 - def `extendLabels` ()
 - Reveives index to slice mask + label id --> returns highlighted mask as PNG to client.*
 - def `getMask` ()
 - Reveives index to slice/mask --> returns mask stored on flask server as PNG to client.*
 - def `highlightOrgan` ()
 - Reveives index of slice + x,y coordinates --> returns highlighted mask as PNG to client.*
 - def `initialize` ()
 - Receive Paths to ordered slices, caches slices.*
 - def `labelOrgan` ()
 - Reveives index to slice mask + label id --> returns highlighted mask as PNG to client.*
 - def `postPickleGetMask` ()
 - Receives path to pickle file --> returns mask as PNG to client.*
 - def `saveMasks` ()
 - Reveives path, saves all stored masks as pickle files to path --> returns output path.*
 - def `updateMask` ()
 - Receives index of slice + slider values --> returns updated mask as PNG to client.*

Variables

- `app` = Flask(__name__)
- int `FLASK_PORT` = 4041
- string `FLASK_HOST` = '0.0.0.0'
- `host`
- dictionary `patients` = {}
 - Dictionary which will hold for each patientID a RadiPopGUI object.*
- `port`

5.3.1 Detailed Description

Flask server for RadiPOP segmenter: Bridge between frontend and backend.

5.3.2 Description

Requests to the flask server follow the HTTP protocol. All requests should contain a JSON formatted message as content. The reply from the Flask server is also always a JSON formatted string. The functions/routes of the flask server don't take any arguments in the classical sense, but expect their "arguments" to be provided in JSON format by the client's POST request. Since there is not better alternative, they appear as the function's parameter in the Doxygen documentation.

5.3.3 Author(s)

- Created by Lorenz Perschy (2021, 2022)

5.3.4 Function Documentation

5.3.4.1 correctPartition()

```
def segmenter_flask_API.correctPartition ( )
```

Reveives index to slice/mask + coordinates--> returns partion corrected mask as PNG to client.

Parameters

<i>patientID</i>	The ID of the patient
<i>index</i>	The index of the slice for which to mask should be updated
<i>coordinates</i>	array of coordinates of the form [x0,y0,x1,y1,...,xn,yn]

Returns

JSON: {mask: byte stream} mask as transparent PNG as byte stream

Note: The coordinates array will be used to generate a line that cuts/divides the segmented organs.

Example handling of return image stream in js:

```
fetch(RadiPOP_states.FLASK_SERVER+"/labelOrgan", {
  method: 'POST',
  headers: { 'Content-Type': 'application/json' },
  body: JSON.stringify(data)
})
.then(function(response){ return response.json();})
.then(function(data){
  bytestring = data["mask"];
  img = bytestring.split('\')[1]
```

5.3.4.2 dcm2png()

```
def segmenter_flask_API.dcm2png ( )
```

Receive Paths to dcm files, converts them to PNG.

Included metadata IDs: "PatientID","PatientBirthDate","PatientName", "PatientAge","PatientSex","PatientName","↵
SliceThickness","StudyID","ContentDate"

Parameters

<i>paths</i>	An array with the paths to the slices
<i>low_clip</i>	lowest pixel value (Recommended:850)
<i>high_clip</i>	highest pixel value (Recommended: 1250)

Returns

JSON: {message: message, metadata: dictionary}

5.3.4.3 dcm2pngPreview()

```
def segmenter_flask_API.dcm2pngPreview ( )
```

Receive Paths to dcm file, converts it to PNG.

Included metadata IDs: "PatientID", "PatientBirthDate", "PatientName", "PatientAge", "PatientSex", "PatientName", "← SliceThickness", "StudyID", "ContentDate"

Parameters

<i>path</i>	The paths to the dcm file
<i>low_clip</i>	lowest pixel value (Recommended:850)
<i>high_clip</i>	highest pixel value (Recommended: 1250)

Returns

JSON: {slice: slice stream} slice as PNG as byte stream

5.3.4.4 drawOnMask()

```
def segmenter_flask_API.drawOnMask ( )
```

Reveives index to slice/mask + x,y coordinates --> returns drawn on mask as PNG to client.

Parameters

<i>patientID</i>	The ID of the patient
<i>index</i>	The index of the slice for which to mask should be updated
<i>coordinates</i>	array of coordinates of the form [x0,y0,x1,y1,...,xn,yn]

Returns

JSON: {mask: byte stream} mask as transparent PNG as byte stream

Note: The coordinates array will be used to draw a line on the mask.

Example handling of return image stream in js:

```
fetch(RadiPOP_states.FLASK_SERVER+"/labelOrgan", {
  method: 'POST',
```

```

    headers: { 'Content-Type': 'application/json' },
    body: JSON.stringify(data)
  })
  .then(function(response) { return response.json(); })
  .then(function(data) {
    bytestring = data["mask"];
    img = bytestring.split('\\')[1]
  })

```

5.3.4.5 extendLabels()

```
def segmenter_flask_API.extendLabels ( )
```

Reveives index to slice mask + label id --> returns highlighted mask as PNG to client.

Parameters

<i>patientID</i>	The ID of the patient
<i>index</i>	The index of the slice for which to mask should be updated
<i>left</i>	Extend labeling up to index-left
<i>right</i>	Extend labeling up to index+right

Returns

JSON: {left_most_idx: idx, right_most_idx: idx}

Note: The *left_most_idx* and *right_most_idx* correspond to the indices of the slices up to which the labeling has been extended. After the the function has finished use the function API's function `/getMask` to update the masks in your GUI. Example in js:

```

for (let index=parseInt(data["left_most_idx"]); index<parseInt(data["right_most_idx"])+1; index++) {
  $.post(FLASK_SERVER+"/getMask", {
    javascript_data: JSON.stringify({patienID: id, index: idx})
  })
}

```

5.3.4.6 getMask()

```
def segmenter_flask_API.getMask ( )
```

Reveives index to slice/mask --> returns mask stored on flask server as PNG to client.

Parameters

<i>patientID</i>	The ID of the patient
<i>index</i>	The index of the slice for which to mask should be updated

Returns

JSON: {mask: byte stream} mask as transparent PNG as byte stream

Example handling of return image stream in js:

```
fetch(RadiPOP_states.FLASK_SERVER+"/labelOrgan", {
  method: 'POST',
  headers: { 'Content-Type': 'application/json' },
  body: JSON.stringify(data)
})
.then(function(response){ return response.json();})
.then(function(data){
  bytestring = data["mask"];
  img = bytestring.split('\')[1]
```

5.3.4.7 highlightOrgan()

```
def segmenter_flask_API.highlightOrgan ( )
```

Reveives index of slice + x,y coordinates --> returns highlighted mask as PNG to client.

Parameters

<i>patientID</i>	The ID of the patient
<i>index</i>	The index of the slice for which to mask should be updated
<i>x</i>	relative x coordinates (0<=x<=1)
<i>y</i>	relative y coordinates (0<=y<=1)

Returns

JSON: {mask: byte stream} mask as transparent PNG as byte stream

Example handling of return image stream in js:

```
fetch(RadiPOP_states.FLASK_SERVER+"/labelOrgan", {
  method: 'POST',
  headers: { 'Content-Type': 'application/json' },
  body: JSON.stringify(data)
})
.then(function(response){ return response.json();})
.then(function(data){
  bytestring = data["mask"];
  img = bytestring.split('\')[1]
```

5.3.4.8 initialize()

```
def segmenter_flask_API.initialize ( )
```

Receive Paths to ordered slices, caches slices.

Parameters

<i>patientID</i>	The ID of the patient. Must be unique!
<i>paths</i>	An array with the paths to the slices

Returns

JSON: {message: message}

Note: Paths to slices !!!MUST BE ORDERED!!! 0,1,...,n

5.3.4.9 labelOrgan()

```
def segmenter_flask_API.labelOrgan ( )
```

Reveives index to slice mask + label id --> returns highlighted mask as PNG to client.

Parameters

<i>patientID</i>	The ID of the patient
<i>index</i>	The index of the slice for which to mask should be updated
<i>label</i>	Label of organ (1 for liver, 2 for spleen, 0 nothing, >2 other organ)

Returns

JSON: {mask: byte stream} mask as transparent PNG as byte stream

Example handling of return image stream in js:

```
fetch(RadiPOP_states.FLASK_SERVER+"/labelOrgan", {
  method: 'POST',
  headers: { 'Content-Type': 'application/json' },
  body: JSON.stringify(data)
})
.then(function(response){ return response.json();})
.then(function(data){
  bytestring = data["mask"];
  img = bytestring.split('\')[1]
```

5.3.4.10 postPickleGetMask()

```
def segmenter_flask_API.postPickleGetMask ( )
```

Receives path to pickle file --> returns mask as PNG to client.

Parameters

<i>patientID</i>	The ID of the patient
<i>index</i>	The index of the slice the mask refers to
<i>path</i>	The path to the mask file

Returns

JSON: {mask: byte stream} mask as transparent PNG as byte stream

Example handling of return image stream in js:

```
fetch(RadiPOP_states.FLASK_SERVER+"/labelOrgan", {
  method: 'POST',
  headers: { 'Content-Type': 'application/json' },
  body: JSON.stringify(data)
})
.then(function(response){ return response.json();})
.then(function(data){
  bytestring = data["mask"];
  img = bytestring.split('\')[1]
```

5.3.4.11 saveMasks()

```
def segmenter_flask_API.saveMasks ( )
```

Reveives path, saves all stored masks as pickle files to path --> returns output path.

Parameters

<i>patientID</i>	The ID of the patient
<i>index</i>	The index of the slice for which to mask should be updated

Returns

JSON: {outdir: path} path/directory to which the pickle files were written

5.3.4.12 updateMask()

```
def segmenter_flask_API.updateMask ( )
```

Receives index of slice + slider values --> returns updated mask as PNG to client.

Parameters

<i>patientID</i>	The ID of the patient
<i>index</i>	The index of the slice for which to mask should be updated
<i>liver-intensity-slider</i>	Slider value for liver intesity
<i>bone-intensity-slider</i>	Slider value for bone intesity
<i>blood-vessel-intensity-slider</i>	Slider value for blood-vessel intesity

Returns

JSON: {mask: byte stream} mask as transparent PNG as byte stream

Example handling of return image stream in js:

```
fetch(RadiPOP_states.FLASK_SERVER+"/labelOrgan", {
  method: 'POST',
  headers: { 'Content-Type': 'application/json' },
  body: JSON.stringify(data)
})
.then(function(response){ return response.json();})
.then(function(data){
  bytestring = data["mask"];
  img = bytestring.split('\')[1]
```

5.3.5 Variable Documentation

5.3.5.1 app

```
app = Flask(__name__)
```

5.3.5.2 FLASK_PORT

```
int FLASK_PORT = 4041
```

5.3.5.3 FLAST_HOST

```
FLAST_HOST = '0.0.0.0'
```

5.3.5.4 host

```
host
```

5.3.5.5 patients

```
dictionary patients = {}
```

Dictionary which will hold for each patientID a RadiPopGUI object.

Patients are added by the API's /initialize function

5.3.5.6 port

```
port
```

Chapter 6

Class Documentation

6.1 RadiPopGUI Class Reference

Bridge between the flask server/API and the RadiPOP segmenter:

Public Member Functions

- def `__init__` (self, `patient_id`)
Class constructor:
- def `extend_labels` (self, `cur_idx`, `left_extend`, `right_extend`)
Extend labels from current slice to neighbouring slices.
- def `highlight_np_label_array_to_png` (self, `mask_np_array`, `highlight`)
Takes an numpy label array $\text{dim}(n,m,1)$ and returns a RGBA pillow image $\text{dim}(n,m,4)$
- def `highlightOrgan` (self, `slice_idx`, `x`, `y`)
Highlights regions of the mask (organs) that were clicked on by user.
- def `labelMask` (self, `slice_idx`, `newlabel`)
Labels mask at given index at previously selected region.
- def `save_masks` (self, `path`)
Saves masks as pickle file to given path.
- def `slice_dim` (self, `index`)
Returns dimensions of slice images (x,y)

Static Public Member Functions

- def `clip_dcm` (`dcm_file`, `clip_low`=850, `clip_high`=1250)
Read dicom image (.dcm), clips it and returns it as a grey scale PNG.
- def `color_mask` (`mask_np_array`)
Takes an numpy label array $\text{dim}(n,m,1)$ and returns a color mask $\text{dim}(n,m,4)$
- def `correct_partition` (`image`)
Convert PNG mask to 1 channelled label mask.
- def `create_image_stream` (`img`)
Returns base64 bytestream for given input image.
- def `draw_on_image` (`coordinates`, `img`, `correctionMode`=False)
Draws a point or lines on given PNG image.

- def [extract_metadata_from_dcm](#) (dcm_file)
Read dicom image (.dcm) and extract metadata Extracted metadata IDs: "PatientID", "PatientBirthDate", "PatientName", "PatientAge", "PatientSex", "PatientName", "SliceThickness", "StudyID", "ContentDate".
- def [find_organs](#) (img, bones_thresh, blood_vessels_thresh, liver_thresh)
Uses three threshold values to find organs.
- def [np_label_array_to_png](#) (mask_np_array)
Takes an numpy label array dim(n,m,1) and returns a RGBA pillow image dim(n,m,4)
- def [read_pickle_mask_to_np_label_array](#) (path)
Opens mask pickle file and returns it as a np.array.
- def [readPNG](#) (path)
Reads an image (e.g.
- def [update_mask_upon_slider_change](#) (image, bone_intensity, blood_vessel_intensity, liver_intensity)
Sets threshold for liver intensity.
- def [writePillow2PNG](#) (img, outfile)

Public Attributes

- [last_clicked_x](#)
The x-coordinate of the region that was last selected/clicked on.
- [last_clicked_y](#)
The y-coordinate of the region that was last selected/clicked on.
- [masks](#)
Dictionary: key: mask index, value: mask numpy array dim(n,m,1)
- [pathToSlices](#)
List containing the path to the png slice files
- [patient_id](#)
ID of the patient.
- [selected_pixel_value_of_label_mask](#)
The label of the region that was last selected/clicked on.
- [sliceCache](#)
Dictionary: key: slice index, value: slice PNG.

Static Public Attributes

- list [LABELS](#) = [[SPLEEN_LABEL](#), [SPLEEN_LABEL](#)]
list with labels
- int [LIVER_LABEL](#) = 1
Regions in mask with this label value are considered liver.
- int [SPLEEN_LABEL](#) = 2
Regions in mask with this label value are considered spleen.

6.1.1 Detailed Description

Bridge between the flask server/API and the RadiPOP segmenter:

Note:

- For each patient one object of this class should be instantiated.
- An object of this class contains all the slices and masks associated with the patient
- This class also contains static utility functions
- This class is the bridge between the flask server/API and the RadiPOP segmenter (segmentations_utils)

6.1.2 Constructor & Destructor Documentation

6.1.2.1 `__init__()`

```
def __init__ (
    self,
    patient_id )
```

Class constructor:

Parameters

<i>patient_id</i>	The ID of the patient. Must be unique!
-------------------	----------------------------------------

Note:

- For each patient one object of this class should be instantiated.
- An object of this class contains all the slices and masks associated with the patient

6.1.3 Member Function Documentation

6.1.3.1 `clip_dcm()`

```
def clip_dcm (
    dcm_file,
    clip_low = 850,
    clip_high = 1250 ) [static]
```

Read dicom image (.dcm), clips it and returns it as a grey scale PNG.

```
@param dcm_file Path to .dcm file
@param clip_low lowest pixel value
@param clip_high highest pixel value

@return tuple(L (grey scale) Pillow Image, slice index)
```

6.1.3.2 `color_mask()`

```
def color_mask (
    mask_np_array ) [static]
```

Takes an numpy label array dim(n,m,1) and returns a color mask dim(n,m,4)

Parameters

<i>mask_np_array</i>	numpy label array dim(n,m,1)
----------------------	------------------------------

Returns

color mask image dim(n,m,4) --> RGBA

- Default liver color is red (LIVER_LABEL)
- Default spleen color is blue (SPLEEN_LABEL)
- Default for other regions is green

6.1.3.3 correct_partition()

```
def correct_partition (  
    image ) [static]
```

Convert PNG mask to 1 channelled label mask.

Parameters

<i>image</i>	Image (e.g.: RGBA PNG pillow)
--------------	-------------------------------

Returns

label mask

6.1.3.4 create_image_stream()

```
def create_image_stream (  
    img ) [static]
```

Returns base64 bytestream for given input image.

Parameters

<i>img</i>	Image (e.g. Pillow PNG)
------------	-------------------------

Returns

img_base64 stream

6.1.3.5 draw_on_image()

```
def draw_on_image (
    coordinates,
    img,
    correctionMode = False ) [static]
```

Draws a point or lines on given PNG image.

Parameters

<i>coordinates</i>	List of the form [x0,y0,x1,y1,...,xn,yn]
<i>img</i>	Image (e.g.: RGBA PNG pillow)
<i>correctionMode</i>	True/False If true the drawn line acts as an eraser (dividing organs). If false a colored line is drawn on the image. DEFAULT: False

The modifications are made directly on the provided image. No return value

6.1.3.6 extend_labels()

```
def extend_labels (
    self,
    cur_idx,
    left_extend,
    right_extend )
```

Extend labels from current slice to neighbouring slices.

Parameters

<i>cur_idx</i>	Index of reference slice
<i>left_extend</i>	Number of slices to extend the labeling to below reference slice
<i>right_extend</i>	Number of slices to extend the labeling to above reference slice

Returns

(left_most_idx,right_most_idx) The index of the outermost slices to which the labeling was extended

Extends labels left and right from current slice How far the labels are extended is taken from left and right expansion bounds

6.1.3.7 extract_metadata_from_dcm()

```
def extract_metadata_from_dcm (
    dcm_file ) [static]
```

Read dicom image (.dcm) and extract metadata Extracted metadata IDs: "PatientID","PatientBirthDate","PatientName", "PatientAge","PatientSex","PatientName","SliceThickness","StudyID","ContentDate".

Parameters

<i>dcm_file</i>	Path to .dcm file
-----------------	-------------------

Returns

dictionary with metadata information

6.1.3.8 find_organs()

```
def find_organs (
    img,
    bones_thresh,
    blood_vessels_thresh,
    liver_thresh ) [static]
```

Uses three threshold values to find organs.

Parameters

<i>image</i>	Image (e.g.: RGBA PNG pillow)
<i>bones_thresh</i>	bones threshold: [threshold, square_size , min_size]
<i>blood_vessels_thresh</i>	blood vessels threshold: [threshold, square_size , min_size]
<i>liver_thresh</i>	liver threshold: [threshold, square_size , min_size]

Returns

New labelled mask (same size as slice)

The algorithm is:

- After some smoothing, remove every pixel above bones threshold from the image.
- After some smoothing, remove every pixel above blood vessel threshold.
- Everything that then remains above liver threshold is called an organ.
- Use contiguous area divisions to roughly split into organs.

6.1.3.9 highlight_np_label_array_to_png()

```
def highlight_np_label_array_to_png (
    self,
    mask_np_array,
    highlight )
```

Takes an numpy label array dim(n,m,1) and returns a RGBA pillow image dim(n,m,4)

Parameters

<i>mask_np_array</i>	numpy label array dim(n,m,1) @highlight Highlight regions where highlight==label in brighter color
----------------------	----------------------------------------------------------------------------------------------------

Returns

RGBA pillow image dim(n,m,4)

Turns a labelled mask into a transparent (RGBA) PNG.

- Default liver color is red (LIVER_LABEL)
- Default spleen color is blue (SPLEEN_LABEL)
- Default for other regions is green

6.1.3.10 highlightOrgan()

```
def highlightOrgan (
    self,
    slice_idx,
    x,
    y )
```

Highlights regions of the mask (organs) that were clicked on by user.

Parameters

<i>slice_idx</i>	Index of mask/slice to be highlighted
<i>x</i>	x-coordinate of slice (in pixels)
<i>y</i>	y-coordinate of slice (in pixels)

Returns

tuple(Mask where the region specified by x and y is highlighted in brigther colors, pixel_value)

6.1.3.11 labelMask()

```
def labelMask (
    self,
    slice_idx,
    newlabel )
```

Labels mask at given index at previously selected region.

Parameters

<i>slice_idx</i>	Index of mask/slice to be labelled
<i>label</i>	Label to be assigned to previously selected region (either LIVER_LABEL or SPLEEN_LABEL or -1 (for removing label))

Returns

Mask with new label

Note: It is expected that the client has before highlighted an organ with the function [self.highlightOrgan\(\)](#). This determines the region/organ that will be labelled with label.

6.1.3.12 np_label_array_to_png()

```
def np_label_array_to_png (
    mask_np_array ) [static]
```

Takes an numpy label array dim(n,m,1) and returns a RGBA pillow image dim(n,m,4)

Parameters

<i>mask_np_array</i>	numpy label array dim(n,m,1) @highlight Highlight regions where highlight==label in brighter color
----------------------	----------------------------------------------------------------------------------------------------

Returns

RGBA pillow image dim(n,m,4)

Turns a labelled mask into a transparent (RGBA) PNG.

- Default liver color is red (LIVER_LABEL)
- Default spleen color is blue (SPLEEN_LABEL)
- Default for other regions is green

6.1.3.13 read_pickle_mask_to_np_label_array()

```
def read_pickle_mask_to_np_label_array (
    path ) [static]
```

Opens mask pickle file and returns it as a np.array.

Parameters

<i>path</i>	Path to pickle file
-------------	---------------------

Returns

numpy array of mask

6.1.3.14 readPNG()

```
def readPNG (
    path ) [static]
```

Reads an image (e.g.

: PNG file) to numpy array

Parameters

<i>path</i>	Path to image
-------------	---------------

Returns

numpy array of image

6.1.3.15 save_masks()

```
def save_masks (
    self,
    path )
```

Saves masks as pickle file to given path.

Parameters

<i>path</i>	Path to which mask files should be written
-------------	--------------------------------------------

Masks are written as .p (pickle) files

6.1.3.16 slice_dim()

```
def slice_dim (
    self,
    index )
```

Returns dimensions of slice images (x,y)

Returns

(x,y) Dimensions of slices

6.1.3.17 update_mask_upon_slider_change()

```
def update_mask_upon_slider_change (
    image,
    bone_intensity,
    blood_vessel_intensity,
    liver_intensity ) [static]
```

Sets threshold for liver intensity.

Parameters

<i>image</i>	Image (e.g.: RGBA PNG pillow)
<i>bones_thresh</i>	bones threshold
<i>blood_vessels_thresh</i>	blood vessels threshold
<i>liver_thresh</i>	liver threshold

Returns

New labelled mask

Steps:

- Sets thresholds for current slice
- Runs self.find_organs on current slice with new thresholds

6.1.3.18 writePillow2PNG()

```
def writePillow2PNG (
    img,
    outfile ) [static]
```

6.1.4 Member Data Documentation

6.1.4.1 LABELS

```
list LABELS = [SPLEEN_LABEL, SPLEEN_LABEL] [static]
```

list with labels

6.1.4.2 last_clicked_x

```
last_clicked_x
```

The x-coordinate of the region that was last selected/clicked on.

6.1.4.3 last_clicked_y

```
last_clicked_y
```

The y-coordinate of the region that was last selected/clicked on.

6.1.4.4 LIVER_LABEL

```
int LIVER_LABEL = 1 [static]
```

Regions in mask with this label value are considered liver.

6.1.4.5 masks

```
masks
```

Dictionary: key: mask index, value: mask numpy array dim(n,m,1)

6.1.4.6 pathToSlices

```
pathToSlices
```

List containing the path to the png slice files

6.1.4.7 patient_id

```
patient_id
```

ID of the patient.

6.1.4.8 selected_pixel_value_of_label_mask

`selected_pixel_value_of_label_mask`

The label of the region that was last selected/clicked on.

6.1.4.9 sliceCache

`sliceCache`

Dictionary: key: slice index, value: slice PNG.

6.1.4.10 SPLEEN_LABEL

`int SPLEEN_LABEL = 2 [static]`

Regions in mask with this label value are considered spleen.

The documentation for this class was generated from the following file:

- `/Users/lorenz/Desktop/radipop/electron-react/segmenter_flask_API/utility/radipop_gui.py`

Chapter 7

File Documentation

7.1 /Users/lorenz/Desktop/radipop/electron-react/segmenter_flask_API/segmenter_flask_API.py File Reference

Namespaces

- namespace `segmenter_flask_API`
Flask server for RadiPOP segmenter: Bridge between frontend and backend.

Functions

- def `correctPartition ()`
Reveives index to slice/mask + coordinates --> returns partion corrected mask as PNG to client.
- def `dcm2png ()`
Receive Paths to dcm files, converts them to PNG.
- def `dcm2pngPreview ()`
Receive Paths to dcm file, converts it to PNG.
- def `drawOnMask ()`
Reveives index to slice/mask + x,y coordinates --> returns drawn on mask as PNG to client.
- def `extendLabels ()`
Reveives index to slice mask + label id --> returns highlighted mask as PNG to client.
- def `getMask ()`
Reveives index to slice/mask --> returns mask stored on flask server as PNG to client.
- def `highlightOrgan ()`
Reveives index of slice + x,y coordinates --> returns highlighted mask as PNG to client.
- def `initialize ()`
Receive Paths to ordered slices, caches slices.
- def `labelOrgan ()`
Reveives index to slice mask + label id --> returns highlighted mask as PNG to client.
- def `postPickleGetMask ()`
Receives path to pickle file --> returns mask as PNG to client.
- def `saveMasks ()`
Reveives path, saves all stored masks as pickle files to path --> returns output path.
- def `updateMask ()`
Receives index of slice + slider values --> returns updated mask as PNG to client.

Variables

- `app` = `Flask(__name__)`
- `int FLASK_PORT` = 4041
- `string FLASK_HOST` = '0.0.0.0'
- `host`
- `dictionary patients` = {}
Dictionary which will hold for each patientID a RadiPopGUI object.
- `port`

7.2 /Users/lorenz/Desktop/radipop/electron-react/segmenter_flask_↔ API/utility/radipop_gui.py File Reference

Classes

- class `RadiPopGUI`
Bridge between the flask server/API and the RadiPOP segmenter:

Namespaces

- namespace `radipop_gui`
Contains python class for storing patient information, handling flask requests and calculations.

7.3 /Users/lorenz/Desktop/radipop/electron-react/segmenter_flask_↔ API/utility/segmentation_utils.py File Reference

Namespaces

- namespace `segmentation_utils`
Contains code from the original segmenter: Segmentation algorithm, applying thresholds ...

Functions

- `def add_sobel_edges` (mask, img)
Smooth edges Steps:
- `def draw_region_outlines` (mask)
Color the mask light green.
- `def guess_bounds` (regions_map, reference_map)
Guess the bounds/labels of the region based on reference region Guess the bounds/labels of the region based on reference region (generally neighboring slice).
- `def partition_at_threshold` (img, thresh, square_size, min_size, title=None, show_plot=True)
After some smoothing, calculate new mask for img Steps:
- `def save_partition` (mask, path)

Index

`/Users/lorenz/Desktop/radipop/electron-react/segmenter_flask_API/segmenter_flask_API.py`,
33
`/Users/lorenz/Desktop/radipop/electron-react/segmenter_flask_API/radipop_gui.py`,
34
`/Users/lorenz/Desktop/radipop/electron-react/segmenter_flask_API/unitysegmentation_utils.py`,
34
__init__
 RadiPopGUI, 23
add_sobel_edges
 segmentation_utils, 10
app
 segmenter_flask_API, 20
clip_dcm
 RadiPopGUI, 23
color_mask
 RadiPopGUI, 23
correct_partition
 RadiPopGUI, 24
correctPartition
 segmenter_flask_API, 14
create_image_stream
 RadiPopGUI, 24
dcm2png
 segmenter_flask_API, 14
dcm2pngPreview
 segmenter_flask_API, 15
draw_on_image
 RadiPopGUI, 24
draw_region_outlines
 segmentation_utils, 11
drawOnMask
 segmenter_flask_API, 15
extend_labels
 RadiPopGUI, 25
extendLabels
 segmenter_flask_API, 16
extract_metadata_from_dcm
 RadiPopGUI, 25
find_organs
 RadiPopGUI, 26
FLASK_PORT
 segmenter_flask_API, 20
FLASK_HOST
 segmenter_flask_API, 20
getMask
 segmenter_flask_API, 17
 guess_bounds
 highlight_np_label_array_to_png
 highlightOrgan
 RadiPopGUI, 26
 RadiPopGUI, 27
 segmenter_flask_API, 17
host
 segmenter_flask_API, 20
initialize
 segmenter_flask_API, 17
labelMask
 RadiPopGUI, 27
labelOrgan
 segmenter_flask_API, 18
LABELS
 RadiPopGUI, 30
last_clicked_x
 RadiPopGUI, 30
last_clicked_y
 RadiPopGUI, 31
LIVER_LABEL
 RadiPopGUI, 31
masks
 RadiPopGUI, 31
np_label_array_to_png
 RadiPopGUI, 28
partition_at_threshold
 segmentation_utils, 11
pathToSlices
 RadiPopGUI, 31
patient_id
 RadiPopGUI, 31
patients
 segmenter_flask_API, 20
port
 segmenter_flask_API, 20
postPickleGetMask
 segmenter_flask_API, 18
radipop_gui, 9
RadiPopGUI, 21
 __init__, 23
 clip_dcm, 23

- color_mask, 23
- correct_partition, 24
- create_image_stream, 24
- draw_on_image, 24
- extend_labels, 25
- extract_metadata_from_dcm, 25
- find_organ, 26
- highlight_np_label_array_to_png, 26
- highlightOrgan, 27
- labelMask, 27
- LABELS, 30
- last_clicked_x, 30
- last_clicked_y, 31
- LIVER_LABEL, 31
- masks, 31
- np_label_array_to_png, 28
- pathToSlices, 31
- patient_id, 31
- read_pickle_mask_to_np_label_array, 28
- readPNG, 29
- save_masks, 29
- selected_pixel_value_of_label_mask, 31
- slice_dim, 29
- sliceCache, 32
- SPLEEN_LABEL, 32
- update_mask_upon_slider_change, 29
- writePillow2PNG, 30
- read_pickle_mask_to_np_label_array
 - RadiPopGUI, 28
- readPNG
 - RadiPopGUI, 29
- save_masks
 - RadiPopGUI, 29
- save_partition
 - segmentation_utils, 12
- saveMasks
 - segmenter_flask_API, 19
- segmentation_utils, 10
 - add_sobel_edges, 10
 - draw_region_outlines, 11
 - guess_bounds, 11
 - partition_at_threshold, 11
 - save_partition, 12
- segmenter_flask_API, 12
 - app, 20
 - correctPartition, 14
 - dcm2png, 14
 - dcm2pngPreview, 15
 - drawOnMask, 15
 - extendLabels, 16
 - FLASK_PORT, 20
 - FLASK_HOST, 20
 - getMask, 16
 - highlightOrgan, 17
 - host, 20
 - initialize, 17
 - labelOrgan, 18
 - patients, 20
 - port, 20
 - postPickleGetMask, 18
 - saveMasks, 19
 - updateMask, 19
 - selected_pixel_value_of_label_mask
 - RadiPopGUI, 31
 - slice_dim
 - RadiPopGUI, 29
 - sliceCache
 - RadiPopGUI, 32
 - SPLEEN_LABEL
 - RadiPopGUI, 32
 - update_mask_upon_slider_change
 - RadiPopGUI, 29
 - updateMask
 - segmenter_flask_API, 19
 - writePillow2PNG
 - RadiPopGUI, 30